

All Pupil II

Now Available with LED Technology

All Pupil II is based on the design of the number one selling All Pupil indirect. Look out for the new features designed to make your examinations easier and more comfortable.

Go Wireless (Patented)

That's right, no more wires and cords! The freedom to move from room to room without any constraints. Keeler is delighted to introduce this great "wow" feature to the number one selling indirect ophthalmoscope.

Outstanding Brightness for Difficult Examinations

With 2000 lux of light, which can be turned down to 2% of total illumination when necessary, All Pupil II offers maximum illumination with total control.



P D Range: 47 - 75mm

All Pupil II LED Module

The All Pupil II LED plugs simply into the back of your All Pupil II to give you all the benefits of LED illumination.

- LED provides better detail, brighter, whiter illumination
- Freedom from purchasing bulbs
- The LED life expectancy is up to 10,000 hours
- 5 year warranty



Lightweight and Comfortable

The new design weighs just 520 grams and extra padded materials have been incorporated into the headband to make this indirect extremely comfortable for all examinations.

Simple Controls

The All Pupil II has a unique single control to change aperture size and position in one movement to view through all pupil sizes.

IR and UV barrier

There is a built in full time IR and UV barrier which guarantees improved patient protection and safety.



Longest life, brightest, better detail



25% brighter than standard bulbs



New Slimline Lithium Polymer battery weighs less than 2 ounces



Combined LED and Lithium Polymer battery can provide up to 6 hours of use on a single charge