



iPac™

Handheld Pachymeter

Intuitive, Inspired Pachymetry.

iPac™ Handheld Pachymeter

The Reichert® iPac™ is the most advanced handheld pachymeter available.

This easy-to-use, lightweight, hand held instrument includes all of the features that you have been searching for in other hand held pachymeters, but have been disappointed with your options. The iPac will impress you with features you have only dreamed about... Bluetooth® wireless connectivity, one-button navigation, rotating color LCD screen and rechargeable lithium ion battery.

Only the iPac provides the most accurate CCT measurements by identifying a greater variety of patient statistics during the measurement process. Advanced Signal Processing converts patient data into custom algorithms that yield only the most accurate measurement.

The intuitive user interface, accessed through a single 5-way button, makes setup easy. You can customize the color display to show IOP correction, standard deviation and L/R or OD/OS notation in seconds. With your pre-programmed user settings, exams become more efficient for you and more comfortable for your patients.



iPac Features:

- **Bluetooth Wireless Connectivity** – Send measurement data directly to your EMR system or Bluetooth compatible printer.
- **One Button Control** – All functions are easily accessed through a single, 5-way control button.
- **Rotating Color LCD** – Large, easy to read display automatically rotates to any viewing orientation.
- **Advanced Battery Technology** – Rechargeable lithium ion battery keeps the iPac ready to use. Power-up with the USB cable (included) or Charging Cradle (optional).
- **Intuitive User Interface** – Easy set-up allows you to quickly customize personal performance needs.
- **Ergonomic Design** – Lightweight, balanced design fits comfortably in your hand.
- **3 Year Warranty** – Industry leading 3 year warranty gives you peace of mind.
- **ADA Eligibility** – Purchase of your iPac may qualify you for ADA tax benefits.

Specifications

Measurement Range	200 - 1000 microns
Accuracy	± 5 microns
Resolution	1 micron
Dimensions	6.69 in. L x 1.77 in. W (17 cm L x 4.5 cm W)
Weight	3.53 oz (100 g)
Power	Lithium Ion Battery

Catalog Numbers

16040	iPac Pachymeter – <i>Includes Carrying Case and USB Wall Charger</i>
16041	iPac Charging Cradle
16042	iPac Rechargeable Lithium Ion Battery
16043	Portable Bluetooth Printer

Proudly Made in the USA.

Visit Reichert Technologies at www.reichert.com or contact your Authorized Reichert Technologies Distributor.

Important Notice

For your own protection, only purchase your Reichert Technologies instruments from and have them serviced by an Authorized Reichert Distributor.

16040-110-Rev-A ©2011 Reichert, Inc. Specifications subject to change without notice. Reichert is a registered trademark of Reichert, Inc. Reichert Technologies, "Advancing Eye Care. Preserving Sight." and iPac are trademarks of Reichert, Inc. Bluetooth word mark and logos are registered trademarks of Bluetooth SIG, Inc.

Reichert
TECHNOLOGIES
Advancing Eye Care. Preserving Sight.™

Reichert, Inc.
3362 Walden Avenue, Depew, NY 14043
Toll Free: 888-849-8955 • Tel: (716) 686-4500
Fax: (716) 686-4545 • www.reichert.com

ISO 9001 Certified
Reichert products are designed and manufactured under quality processes meeting ISO 9001 requirements.