

Professional Ophthalmoscope

Total Control at Your Fingertips

The Keeler Professional Ophthalmoscope has been designed to ensure ultimate control. Lenses, graticules and filters can all be introduced at the touch of a fingertip. No need to pull away from the patient, the Professional Ophthalmoscope makes routine examination far easier.

Xenon Illumination

We recognise that illumination is of fundamental importance to good diagnosis. That's why we have selected Xenon illumination for our ophthalmoscopes. This brighter and whiter longer lasting light will meet any diagnostic challenge! With Xenon you have the power to diagnose even through cloudy media. And because it is longer lasting, that means fewer bulb changes giving greater peace of mind.

Positive Action Maywheel

The positive action easy access maywheel not only allows you to dial up the required lenses quickly, it also allows you to know where you are without the need to pull away from the patient.

Comprehensive Lens Range

The extensive range of +29D to -30D in single dioptre steps ensures easy examination of all ocular structures and covers a wide range of patient/user prescriptions. A quick step swing in +/- 20D helps speed up diagnosis.

Swing Over Filter

The red free filter for enhanced vessel examination can be used in conjunction with any of the 6 diagnostic beams.

Brow Rest

The integral soft brow rest allows for easy positioning and protects your lenses if you are a spectacle wearer.



Part Numbers

- | | |
|-------------|---|
| 1125-P-1180 | 3.6v Professional Ophthalmoscope on Lithium Handle (no stand adaptor) |
| 1125-P-1124 | 3.6v Professional Ophthalmoscope on Lithium Handle with stand adaptor |
| 1128-P-1004 | 3.6v Professional Ophthalmoscope head only |



Specialist Ophthalmoscope

The Choice of the Ophthalmic Profession

The Keeler Specialist offers precision optics, versatile features and practical ergonomics making it the number one choice for ophthalmic practitioners.

Unequaled Quality

The Specialist ophthalmoscope is perfectly balanced for precise ocular examination. The comprehensive lens range (+44D to -45D in 1D steps) ensures clear and accurate diagnosis even in the most difficult eyes. The Specialist lens range is made up of tiny individual ground glass lenses. The fine gearing of the Morton Track lens system allows you to dial up any lens quickly and easily.

Versatility

The graticule range offers a total of six diagnostic beams for comprehensive ocular examinations. In addition to a choice of three viewing beam sizes the Specialist includes a red free filter for vessel examination, a fixation cross for

assessment of eccentric fixation and a slit to determine changes in retinal elevation. A snug orbital fit optimises the field of view and guarantees increased comfort for both practitioner and patient.

Pupillometer

The pupillometer comparison disc allows quick assessment of pupil size.



Part Numbers

- | | |
|-------------|---|
| 1132-P-1025 | 3.6v Specialist Ophthalmoscope on Lithium Handle (no stand adaptor) |
| 1132-P-1025 | 3.6v Specialist Ophthalmoscope on Lithium Handle with stand adaptor |
| 1132-P-1050 | 3.6v Specialist Ophthalmoscope head only |
| 1132-P-5000 | 2.8v Specialist Ophthalmoscope head only |

