



Professional Combi Retinoscope

Streak or Spot – just a Bulb Change!

There are those that love Streak and those that love Spot. The Keeler Professional Combi offers the best of both worlds. The unique optical design provides the first truly bimodal Retinoscope with just a quick bulb change.

Double Neutralization Check

The Professional Combi is the only Retinoscope which combines the ability to switch from Spot to Streak with the facility to double check neutralisation in both modes.

Easy Access Control

The Professional Combi uses rare earth magnet technology to provide optical precision refraction. A single control is used for rotation, convergence, and divergence. A choice of apertures to optimise brightness or resolution of the retinal reflex in 4mm or 1.7mm The Professional Combi comes with easily interchangeable soft brow rest for spectacle wearers.

Light Output

With up to 600 lux of brilliant illumination, the Professional Combi Retinoscope guarantees you more than enough light to refract even the most difficult eyes.

Near Point Cards

Fixation cards for Dynamic Retinoscopy simply clip on and off for easy use.

Part Numbers

1322-P-1060

3.6v Professional Combi Streak Retinoscope on Lithium Handle (no stand adaptor)

1302-P-1007

3.6v Professional Combi Streak Retinoscope head only

Finoff Transilluminator

For use when the presence of intraocular tumors in patients with cloudy media is suspected. The Finoff transilluminator illuminates the fundus via the sclera and any tumors present in the light path will dim the normally bright red fundal reflex as seen through the pupil.

Part Number

2903-P-1126

3.5v Transilluminator (Head only)

2903-P-1118

2.8v Transilluminator (Head only)

