



**3000**

Manual Pneumatic Exam and Treatment Chair

**Tradition and Innovation** – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.



# 3000 Pump Hydraulic Lift and Manual Recline Chair

## Ergonomic Features

- Adjustable padded headrest locks in place with a single hand.
- Foot bar and armrests safely support patient's full weight.
- Armrests rotate to allow patients in or out from both sides of the chair.
- Rotation lock allows chair to swivel 360° for improved patient access.

## Pump Hydraulic Base for Effortless Raise

- Hydraulic mechanism allows you to easily raise the patient by pumping the 3000's foot pedal. Automatically lowers when pedal is moved into the down position.

## Manual Reclining Top

- Chair back reclines up to 90° with a single hand. Convenient release mechanism on both sides of the chair.

## Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Removable cushions for easy cleaning or replacement.
- Interchangeable seat and back cushions.
- Choose from 14 standard vinyl colors and endless custom selections – call Reliance at 1-800-735-0357 for a color card.

## Tech Specs

- Top and base shipped together for easy installation. Base height: 28" raised, 19" lowered. Base dimensions: 23" diameter. Overall width: 25". Depth upright: 33 1/2". Depth reclined: 42". Seat dimensions: 19" by 19". Armrest distance: 20". Shipping weight: 220 lbs.

