

Name: ID: 1111111111 DOB: 14-01-1971




24-2 ZETA Fast

Fixation target: central
 Gaze errors: ---
 HK: 2/9 !
 FPOS: 0/12
 FNEG: 9/11 !!
 Duration: 06:00
 Fovea: ---


Stimulus: Ill, White
 Background: W:31.5 ASB
 Strategy: ZETA Fast
 HoV @10deg: 33 dB (-3.8 dB/10°)

Pupil diam.: 3.2 mm
 Rx used: +2.00 DS +0.00 DC x 000°

Date: 14-01-2020
 Time: 10:03:24
 Age: 49




Total Deviation





Pattern Deviation

GHT:
 Outside Normal Limits
 VQi: 69.8%
 MDh: -10.69 dB
 PSD: 12.45
 CPSD: 1.72
 SFh: 11.2
 MS: 21.60 dB



- <math>< 0.5\%</math>
- ▨ <math>< 1\%</math>
- ▩ <math>< 2\%</math>
- <math>< 5\%</math>
- △ >



Exam comments:



Optopol Technology
 Zabia 42
 42-400 Zawiercie
 Poland
 tel. 0048 32 670 91 73

Operator:
 AppUser:
 Device: PTS 2000:7700355;FW:20180913;SW:3.4.0.21
 Copyright 2020 Optopol Technology
 PTS ver.3.4.0.21